

Disposable Pulse Lavage

Indications

It can be used for wound debridement, soft tissue debridement and cleaning of surgical sites.



EP18-01 (Joint Type)

Atomizing nozzle, atomizing nozzle with brush



EP18-02 (Spine Type)

Atomizing nozzle, lengthened atomizing nozzle



EP18-03 (Traumatic Type)

Atomizing Nozzle, Shower Nozzle



Features

- ※ Ergonomic design for easy operation
- ※ One machine with multiple heads, making reasonable clinical choices according to different wound surfaces and surgical requirements
- ※ The flush-integrated spray head is designed to quickly suck the flush away, preventing the spread of contaminants
- ※ The spray head is equipped with a splash guard, which can effectively prevent flushing splashes and protect medical personnel
- ※ High-performance alkaline battery ensures long-term operation of the main engine

Specifications

Code	Type	Components
EP18-01	Joint Type	Handpiece, atomizing nozzle, atomizing nozzle with brush
EP18-02	Spine Type	Handpiece, atomizing nozzle, lengthened atomizing nozzle
EP18-03	Traumatic Type	Handpiece, atomizing Nozzle, Shower Nozzle

www.zkskmedical.com
 TEL: 0086 186 7820 2081 English line * 0086 137 9218 3575 español
 linea *
 Email: summer@bjzksk.com shan@bjzksk.com